InFocus[™] Dynamic Optical Focusing Systems







2530 Frontier Ave., Boulder, CO 80301 • Tel: (303) 440-4544 • Fax: (303) 440-4144 • E-mail: info@infinity-usa.com • Web: www.infinity-usa.com

Copyright 2014 Infinity Photo-Optical Company. All Rights Reserved.



InFocus[™] Dynamic Optical Focusing Systems

INSTRUCTIONS AND TECHNICAL DISCUSSION

- Imparts DYNAMIC OPTICAL FOCUSING to Almost ANY Infinity-Corrected Microscope
- Corrects Spherical Aberration BETTER than Correction Collared Objectives
- Works with ALL Infinity-Corrected Objectives
- Color-matched for Zeiss, Leica, Nikon or Olympus—and others
- Allows YOU to Test Your Objective's Performance
- Breadboard or Stand Options

InFocus is a unique optical device—in fact, a system—that allows you to realize the full potentials of optical microscopy. InFocus can be thought of as the "ultimate" InfiniTube when infinity-corrected objectives are used with it. But where Infinity Photo-Optical's simpler InfiniTube Series lenses impart minor focus adjustment, InFocus imparts true dynamic optical focusing. Whether you use InFocus as a self-contained microscope in a breadboard setup—or attached to a microscope stand—it provides added degrees of convenience and perfection. True optical focusing means that a microscope equipped with InFocus can be deliberately defocused so that a new position can be set for the relative working distance of microscope object to objective. This is more significant than may at first be thought.

Beyond its simpler industrial and general uses, InFocus is a significant breakthrough wherever cover glasses must be used. In biological and high-tech industrial applications, this means that it can be used to diminish or correct the problem of *preparation induced* spherical aberration (SA). (In SA, an optical system cannot focus all the light at a common plane resulting in a fogging of image details). This is the single most overlooked reason for poor microscope imagery. Even the finest optics can be thrown out of correction when cover glasses, thick or highly-refractive media are used.

This is all due to the fact that the relative position of the objective's focus is determined by the thickness of the cover glass (or lack thereof) and the medium in which the subject is immersed. Manufacturers *claim* to correct their biological or special industrial objectives for use with a *defined* cover glass thickness with the preparation just beneath it. Deviating from this factory setting creates (for various optical reasons) SA whereby image detail is never properly focused or captured. The resulting image can even become useless blurring—particularly in fluorescence applications. But, unless a perfect objective is used *perfectly* with a cover glass (biologically, usually set for 0.17mm) designed for it—with a specimen mounted just beneath in a mounting medium of specific refractive index—SA and a degraded image—are inevitable. This is why most photomicrographs are imperfect. As a consequence, many published images are more or less "fuzzy" compared to how they would be if captured with a spherically corrected system.

Background of InFocus. Until InFocus, there were two ways to correct SA in microscopes. One was to push in or pull out the drawtube of old-fashioned monocular microscopes. For every 10mm of movement, approximately 0.01mm of cover glass correction could be made. This was not very efficient, especially with thicker preparations where only about 20mm of down movement was even possible. The other way was to use microscope objectives fitted with correction collars that separated the lens elements and reset the focus obtained. But this method—still used today—is cumbersome, difficult and usable only with objectives so-equipped to begin with. Objectives which do not incorporate correction collars can never be utilized to advantage. More, just by touching the objective to make the correction, decentering or compression of the sample are almost assured.

The original InFocus eventually went on to be motorized and computer-activated in order to determine its bestcontrast setting, pricing itself into a small, exclusive market. Although simpler systems were offered, their functional ranges needed supplemental lenses to extend operational efficiency where needed. Still, the original InFocus was such an advance, that it went on to win *two* R&D 100 awards.

In recent years, Infinity Photo-Optical introduced the K-Series long-distance microscopes equipped with patented CentriTel® focusing. CentriTel is an internal focusing system which does so without essential magnification change. This allowed us to revisit the InFocus concept and apply features derived from CentriTel to it. During that time, camera sensor technology advanced considerably. As a result, it is now possible to offer an *improved* **InFocus Module that is simpler, more universal—and does not require computer-activation in order to determine its best setting**—putting the price and image enhancements of InFocus within reach of most research and industrial microscopists. In addition, the improved InFocus has minimal magnification change through its focal translation and centers on a c.1x factor.

What InFocus is. InFocus is an integrated microscope *tube with patented optics*, complete with a tube lens for use with infinity-corrected objectives of various makes (more on this later). When mated to an infinity-corrected microscope objective (either integrated in a breadboard setup or on a microscope stand), InFocus is a complete microscope—but with one very considerable difference: it has an internal focusing adjustment which is actually *finer* than any mechanical focus can be. A dynamic optical focusing system. Turning a knurled ring activates InFocus to translate focus from above, through and below the object. For industrial setups, InFocus allows focal translation and "touch up" even if the unit is bolted in place. For biological or high-tech industrial applications, it goes one step further by utilizing its ability to alter tube length—as a drawtube would—but without changing outer dimensions. Also, its operative range is far, far greater than that afforded by any drawtube. As a result, InFocus imparts the ability to correct spherical aberration to *any objective used with it*—even those already equipped with correction collars (which should then be set at their null position).

By adding InFocus' dynamic optical focusing to any microscope—either onto the stand or in a breadboard system—it becomes possible to do things that previously required special equipment. Now, *every* microscope is *micromanipulator-ready*. Now, every microscope can be equipped *beyond the capabilities of what correction collars*—*which do NOT provide dynamic optical focusing*—*can do*. In fact, the degree of correction possible from InFocus is considerably greater than that of typical correction collar-equipped objectives. And, if quality results *cannot* be obtained at a reasonable setting of the InFocus, your objective may be a defective one.

The improved InFocus has an operating range approximately *four times* that of the original. This opens new vistas for use. Now, it is possible to set the initial SA correct position and, from that starting point, use the *optical* focus instead of the mechanical to select various Z-depth planes. There are now two ways to use the InFocus for specific results. The first is to set it for SA correction. For single one-plane image capture, that is sufficient. But for Z-stacking (with software such as Helicon FocusTM or ZareneTM), the SA correction can be made and then reestablished at another depth obtained by using the mechanical focus. That is the traditional way to use InFocus. However, it is now possible—due to the extended range of the improved module—to set SA correction and then use *only the InFocus to dynamically select different depths—abandoning the use of mechanical focus altogether.* This feature of the improved module makes it possible to focus *without moving the objective* through delicate membranes which might otherwise be subject to medium compression. What is centered to start is still centered at the finish of a focal sweep.

Regardless of the method used, a microscope equipped with InFocus can provide images better than one which is not. It is as simple as that. But it was neither simple nor easy to make that *happen* during the development of InFocus.

DYNAMIC FOCUSING: SUPERIOR TO CORRECTION COLLARS

By 1830, both Amici and Lister had independently observed that spherical aberration was induced by using cover glasses with high apertured microscope objectives. This severely diminished image quality. In 1837, Andrew Ross—acting on Lister's suggestion—made the first objective with a built-in correction collar. The correction collar (and some derivatives of it such as Watson's Jackson Tube Length Corrector

which transferred the action to a separate lens system in the finite beam) has remained the only means to compensate *preparation-induced* aberrations ever since.

With correction collars, the objective's own lens system can be actively separated. At a normal null position (with biological objectives, generally 0.17mm cover/medium thickness), a correctly made preparation will provide a "best image." If the cover/medium thickness deviates from the norm, the correction collar is turned to set a new "best image" which is refocused by using the mechanical focus. Turning the correction collar itself does not impart focal *translation*. Consequently, correction collars activate optical changes *within* the objective and *only for that particular objective*. When used with infinity-corrected systems, correction collared objectives operate on the beam characteristics *before they are acted upon by the infinity tube lens system*. The infinity tube lens system remains STATIC.

With InFocus, one common module—usable with virtually every objective—acts almost precisely in reverse. With InFocus, the objective remains STATIC but the infinity tube lens system becomes DYNAMIC. In effect, InFocus "forces" an objective to seek a new frontal focus. Then, the previous plane can be re-established by the mechanical focus with increased fidelity because InFocus has automatically compensated for all factors: the cover glass' thickness, deviations in media refractive indices and refractive errors. Because InFocus imparts DYNAMIC OPTICAL FOCUSING, different depths can be focused upon (once the initial "best image" is obtained). This feature is of immense advantage to those who cannot allow the objective to move in any way, particularly when observing delicate tissues which are easily displaced by compression or when industrial processes can come in contact with the objective.

Since InFocus is a *module* rather than a *feature* built into a single objective, it can be used to critically test the quality of objectives. If the chosen objective is claimed to be factory-set for 0.17mm and the InFocus cannot correct its image, it is probably faulty. And, since the InFocus works with all objectives, you need not buy correction collared objectives if a simpler mounting is available. If all this were not enough, InFocus can be mounted on dovetails so that one, usable with one make of objective can be switched with another InFocus usable with another make. You can actually *choose* which objective to use from among various makes and types—a dream since the early years of microscopy.

Four Major Brands—THREE solutions required. Years ago, it was reasonably simple to use one brand of microscope objective next to another. Makers pledged to use the "RMS Thread" to fit objectives to microscopes and, even if one make used a slightly different tube length, the drawtube could correct that. Most all microscopes were monoculars. Even if one maker's objective needed a different eyepiece to work with it than another's, it was easy to pull one eyepiece out and put another in its place. But with modern binocular microscopes and the advent of flat field (plan) objectives, things got complex. Then, when the advantages of infinity-correction became too obvious to ignore, makers settled upon their own optical solutions. Zeiss chose to keep 160mm tube length and RMS mounting thread (now also offering M27 x 0.75), the same objective length as before and transferred their compensation from "KPL" eyepieces to the infinity tube lens. Leica did similarly, selecting 200mm as the tube length and putting their "Delta" correction in the tube lens. Thus, Zeiss and Leica use "color-specific" tube lenses. Nikon and Olympus chose to use various threads and dimensions (Nikon 200mm and Olympus 180mm tubelength), etc. but achieve all color correction in the objective itself. Nikon and Olympus tube lenses can therefore be termed "neutral." Others also use the "neutral" correction, such as Mitutoyo, Edmund and Infinity Photo-Optical.

In order to offer InFocus as a viable product, it was necessary to replicate the corrections for the two "colorspecific" brands and provide a "neutral" system for all others. This was a major challenge but Infinity succeeded in accomplishing it. Consequently, there are not one—but *three* InFocus modules. The Zeiss and Nikon/Olympus modules are the same overall length while the Leica is 9.5mm longer. The InFocus' focal length has been selected to be c.160mm; therefore, a Zeiss objective operates at its rated magnification, an Olympus at 0.88x its rated magnification, and a Nikon, Mitutoyo or Infinity Photo-Optical objective at 0.80x its rated magnification. However, amplifiers can attach to the InFocus to permit alteration of these "raw" factors. This approach allowed us to set the rear accessory distances to be the same. With the successful development of *three* different InFocus modules—all compatible with a whole line of accessories, the InFocus is properly a *system* unto itself. Bluntly, if your microscope is not equipped with InFocus, it probably is *not in focus*.

QUICK-START CHECKLIST:

NOTE: Be sure to keep the condensing aperture fully-open while setting InFocus. This will allow you to easily see the tremendous increase in contrast and determine the best setting. Modern DSLRs and other digital cameras have provisions to electronically increase screen magnifications to dramatically show this.

- Mount the InFocus via its proper dovetail adapter on a microscope stand or, if used in a breadboard setup, use a T-tube to secure it. Clamps and other accessories permit this to be done (see Drawings). If the bottom part of the InFocus contacts any part of your microscope, a small T-tube should be interfaced—but this is unlikely on microscope stands.
- 2) Common parfocalize the system. You can now absolutely parfocal your microscope. When your microscope is set for the common parfocal position, swinging one objective into place from another results not in a generally focused image, but one which is essentially at the same plane of focus as the other.
- 3) To SET for common parfocality, use a low power (e.g., 10x) objective and then use the highest numerically apertured high dry objective (HDO). Generally, this should be a 40/0.75 or 60/0.90 objective. Focus the 10x and then swing to the HDO. Swing back to the 10x and turn InFocus' adjustment to reset focus. Repeat this. After the 10x has been refocused by the InFocus a second time, switching between the objectives will result in the same plane being in focus, regardless which objective is used. The microscope is now set for common parfocality. The goal is to have your microscope set so that "infinity" is virtually absolute between two objectives whose depth of field v. depth of focus characteristics are near opposite. This establishes the starting point for InFocus' operation.
- 4) Your InFocus was set at the factory so that when the two arrows on its focusing ring/body face each other, it is at infinity. If you wish to reset them, use a 2-56 Allen wrench (2mm metric) to *slightly detent the metal-colored setscrew*, turn, and *gently* lock in place.
- 5) You will note MINUS (-) and PLUS (+) markings on the InFocus body. Minus makes the unit correct for progressively *thinner* preparations and Plus makes the unit correct for progressively *thicker* preparations. Please see the InFocus Graphics for an explanation as to what these mean exactly. In any case, the InFocus can be used as simply as capturing one image, then another, all while moving the InFocus control in several steps and re-establishing the same plane by mechanical means. *InFocus CAN be as simple to use as that: One image will stand out as the best; one image will be the one you wish to "publish."*

Advanced Use. While indiscriminate use of InFocus is possible—as outlined above—InFocus comes into its own when it is used with more technical attention. The improved InFocus now permits *two* techniques to be utilized:

- 1) Mechanical Focus Reset (MFR).
- 2) Dynamic Optical Focus Sweep (DOFS).

Mechanical Focus Reset. Here, the InFocus is progressively turned and a mechanical focus re-established for the plane of interest. The contrast will either improve or deteriorate depending on the match for SA correction. If a thin preparation is used, first start with the Minus side of the translation. If a thick preparation is used, "weight" the translation to the Plus. Observe the increase in contrast at each setting. If you are using InFocus with a

camera that allows electronic magnification increase, this may be very simple to determine. You can focus at different Z-depths and process for stacking or deconvolution, later.

Dynamic Optical Focus Sweep. This technique is possible due to the increased efficiency of the improved InFocus module. The best contrast is found by using MFR. It differs thereafter, since, from *that position, the InFocus is used to seek different Z-depths.* The mechanical focus *is not used at all.*

[The improved InFocus has such a sufficient operating range, that even a diatom mounted in a medium such as Pleurax (1.8 index of refraction) can be brought to greatest contrast at c.90-degrees of the turn to Plus, leaving at least that many degrees of additional turn of the InFocus dial. Even then, if the InFocus is turned still more to the Plus, the image will *degrade*. In effect, the InFocus is then "asking" for a diatom at an even *lower* level to be imaged.]

DOFS is particularly useful where delicate membranes must be imaged and where *any* movement of objective could cause displacement, rendering the collective images captured useless. Since any focal changes from the InFocus are independent of connection or contact with the objective and preparation, it allows images to be captured as never before.

InFocus Graphics. Full details on how the InFocus works, what preparations it can be used with and the SA corrections possible are contained in the InFocus Graphics section attached to the Instructions PDF. The InFocus Graphics also show the relative positions of the objective in different media and/or thickness situations.

Applications and Camera Formats. If you are doing deep cell translations; attempting point-spread functions; 3-D imaging or confocal microscopy; if you are examining wafers and wish to avoid all possible contact; if you are building a breadboard instrument—or simply trying to capture accurate imagery—InFocus will allow you to do things previously thought impossible. **Every microscope now becomes micromanipulator-ready.** And, InFocus has a whole *system* of accessories to back it. With them, it is possible to use virtually *every type of camera that has ever been used in photomicrography.* This holds especially true if you have one of the advanced research-grade C-mount cameras from the major microscope makers. All their software can be utilized.

Modern digital cameras are now so efficient that the improvement in contrast provided by InFocus is easy to see. This is especially true when electronic magnification is increased on the display or to a monitor. The determination of the point where contrast is maximized by InFocus is such, that computer determination (which was once considered a necessity) is no longer required.

The Plus-Weighting C-Mount. Most research-grade cameras offered by the major microscope makers are 2/3in. C-mount types. If deep focal translations, high refractive index media, thick covers, etc., are used or contemplated, the Plus-range of the InFocus becomes more and more necessary to achieve corrections. The Plus-Weighting C-mount can be shortened in two steps of 12mm so that common parfocality occurs with most of the Plus-range still able to be used. If you are planning to do such work, we highly recommend the Plus-Weighting C-mount instead of the regular types we offer.

Care and Cleaning. The InFocus should be treated as the fine optical instrument it is. Care should be taken to keep dust and dirt off external surfaces. If the front element should ever need to be cleaned, use lens tissue moistened by a lens cleaner. NEVER USE SOLVENTS of any kind. The InFocus should NEVER be opened. If contaminants enter any inner parts, consult your dealer or INFINITY about cleaning services.

Warranty Service and Questions. Specific details of the warranty are given on the limited warranty statement. In general, all parts and labor are guaranteed for one full year. Should the InFocus Insert become damaged or need service, return it to your authorized dealer with a letter explaining the problem. If you have additional questions, please do not hesitate to contact your dealer or INFINITY PHOTO-OPTICAL COMPANY directly.

InFocus is a trademark of Infinity Photo-Optical Company. U.S. Patent 7,869,139. Additional patents pending.

InFocus[™] System



TO MAJOR-BRAND STAND/FRONT ACCESSORIES

InFocus[™] PRO System



TO MAJOR-BRAND STAND/FRONT ACCESSORIES

InFocus[™] Objective Options

For Breadboard/Industrial Applications

There is an InFocus Module for most major-brand infinity-corrected objectives:



Trademarks of Infinity Photo-Optical Company

The Plus-Weighting C-mount (PWC)

The PWC allows the InFocus Module to be progressively Plus-Weighted.



InFocus MAIN BODY

Normal Range





Plus-Weighted (one 12mm tube removed)

Fully Plus-Weighted (both 12mm tubes removed)

InFocus[™] Graphics

Infinity Photo-Optical Company Boulder, Colorado Tel 303/4404544 Fax 303/440-4144 WWW.Infinity-USA.com

Conventional Focusing



Situation A: Ideal Conditions. Cover glass = 0.17mm; tube length is correct. Intermediary image is in prescribed position relative to frontal conjugate and media.

Situation B: Erroneous Conditions. Frontal conjugate has changed, as has media relationship. Yet, intermediary image is "forced" to form at same position as before.

Optical Infinity

The common depiction of an infinity beam is incorrect...

... Actually, the conic diameter expands--- though no image forms.

Spherically Correct Cover/Preparation Variations for Infinity-Corrected Objectives





Cover too thick









Immersion



Low index medium



Immersion with high index medium

Corrections were obtained by using InFocus.

Cover/Preparation Variations and Relative Objective Positions

(where the objective would "like to be")



Cover too thick



Cover too thin



High index medium



Low index medium



Immersion

Spherically Incorrect Cover/Preparation Variations for Infinity-Corrected Objectives



High index medium

Low index medium



Immersion

Immersion with high index medium

Arrows show divergence of focusing direction from factory setting.

InFocus Microscope System (What it is)



The InFocus module does not work with the Standard Tube Lens. The InFocus *replaces* it.

Video Formats



Unless otherwise noted all data is given for 1/2" camera sensors. For FOV or magnification in other sensor sizes use the factors listed. For example, a FOV of 25mm listed for a 1/3" camera would be $25 \times 0.75 = 18.75$ mm.

INFINITY PHOTO-OPTICAL COMPANY LIMITED WARRANTY

INFINITY PHOTO-OPTICAL COMPANY hereby warrants its products to be free from defects in workmanship or materials for the warranty period set forth below. INFINITY PHOTO-OPTICAL COMPANY, at its option, shall repair or replace the defective product without cost to the purchaser, and such repair or replacement shall be the full extent of this express limited warranty. INFINITY PHOTO-OPTICAL OPTICAL COMPANY shall not be liable for any other damages either direct or consequential.

This warranty is made to the original purchaser, and is effective only on new equipment purchased from INFINITY PHOTO-OPTICAL COMPANY, or a dealer authorized by INFINITY PHOTO-OPTICAL COMPANY to sell the product.

This warranty is valid only when the product is returned to the authorized dealer from whom it was purchased, or returned directly to INFINITY PHOTO-OPTICAL COMPANY, freight prepaid, with proof of date of purchase.

This warranty does not extend to any defect, malfunction or failure caused by misuse, abuse or the use of the product with equipment for which it may not have been intended. Any unauthorized repair voids this warranty.

The warranty period for all products manufactured by INFINITY PHOTO-OPTICAL COMPANY is five (5) years from date of original purchase. Parts or components made or sourced from other manufacturers shall be solely covered by that manufacturer's warranty.

The warranty contained herein is the only warranty made by INFINITY PHOTO-OPTICAL COMPANY. Any implied warranty of merchantability and/or fitness for a particular purpose is expressly excluded from this warranty. INFINITY PHOTO-OPTICAL COMPANY shall not be liable for any expense, loss, incidental or consequential damages which may arise in connection with the use of this equipment. Recovery under this warranty is limited to repair or replacement of the equipment as provided above.

INFINITY PHOTO-OPTICAL COMPANY reserves the right to modify designs, equipment and accessories without notice.