# K-SERIES LONG-DISTANCE MICROSCOPES: MODEL K1 CentriMax™

Suggested List Prices Effective: April 9, 2018

### K1 CentriMax LONG-DISTANCE MICROSCOPES

The K1 CentriMax is Infinity's premier all-inclusive CentriTel<sup>™</sup> optic. Capable of focus from 3 meters to 60mm, depending on the MX-series objective used, the K1 CentriMax maintains essentially constant magnification when focused from a selected working distance.

The K1 CentriMax can be supplemented with coaxial illumination and infinity-corrected microscope objectives for higher-power uses in the InFocus<sup>™</sup> mode. In the InFocus mode, some slight magnification changes can be expected during internal focus. It then has a 2x factor (for example, a 10x objective then operates as a 20x). The standard S Front can be used with Achrovid<sup>™</sup>, Mitutoyo, Nikon or Olympus infinity-corrected objectives. For use specifically with Zeiss ICS-system color-corrected objectives, the S Front can be exchanged with the Z Front.

For use with sensors larger than 4/3-in., the K1 CentriMax can be supplemented with the NTX 2x amplifier. A camera adapter (C or T2) must be added.

It is easy to order a K1 CentriMax suited to your needs. Simply choose a K1 CentriMax Main Body, Objectives and applicable accessories to configure a complete K1 CentriMax system.

#### K1 CentriMax Main Body

860502	K1 CentriMax Main Body with Variable Iris Diaphragm	\$	4,450.00	
K1 CentriMax ROBUSTO™ Main Body				
865502	K1 CentriMax ROBUSTO™ Main Body with Variable Iris Diaphragm Includes NTX Tube (2x) as part of ROBUSTO™ system	\$	5,200.00	

#### **Camera Mounts**

770507	C-mount Video Adapter	\$	50.00
770509	Orientable/Separable C-Mount Video Adapter	•	425.00
770510	Spare C-mount Top for 770509		195.00
770511	Rotatable/Lockable C-mount		125.00
770568	T2 Adapter – Nikon F mount		12.00
770565	T2 Adapter – Canon EOS		12.00
770566	T2 Adapter - Olympus		12.00
770571	ArriFlex PL mount Adapter		150.00

# **Amplifiers**

690568 220568	DL Tube (2x) (video only) LDL 2x Lockable Power Module Amplifier	\$	195.00 275.00		
MX-Series Objectives					
860170 860171 860172 860173 860174 860175	MX-LR (0.11x-0.58x; 3000-530mm WD) MX-1 (0.56x-0.98x; 590-303mm WD) MX-2 (0.85x-1.28x; 383-233mm WD) MX-3 (1.16x-1.60x; 254-170mm WD) MX-4 (1.42x-1.83x; 223-160mm WD) MX-5 (1.83x-2.13x; 168-127mm WD) MX-6 (CONSTANT 4.27x; 75-60mm WD)	\$	465.00 465.00 525.00 465.00 465.00 465.00		
Achrovid™ Nelsonian™ Objectives (infinity-corrected/f=200mm)					
210301 210302 210303 210304 210305 210306	Achrovid 2x/Nelsonian 0.054/ W.D. 37mm Achrovid 2.5x/Nelsonian 0.064/ W.D. 37mm Achrovid 3.3x/Nelsonian 0.077/ W.D. 37mm Achrovid 4x/Nelsonian 0.09/ W.D. 37mm Achrovid 5x/Nelsonian 0.12/ W.D. 37mm Achrovid 10x/Nelsonian 0.21/ W.D. 17.50mm	\$	410.00 410.00 410.00 410.00 410.00 485.00		
991158 991159 992258 992259 992283 991170 991172 992270 992271 991167	Coaxial Illuminator w/ Collector I (for use with low power objectives). Coaxial Illuminator w/ Collector II (for use with objectives 10x and higher). Basic Coaxial Illuminator w/ LED Flashlight Basic Coaxial Illuminator Right Angle w/ LED Flashlight Right Angle Interface for Coaxial Illuminator. Converts Standard into RA Type. Collector I (to convert units that have Collector II) Collector II (to convert units that have Collector I) Collector for LED Flashlight (extra) LED Flashlight for Coaxial Illuminator. Polarizing Filter Set and Retainers	\$	895.00 895.00 960.00 1,475.00 555.00 235.00 235.00 235.00 25.00 210.00		

## **General Accessories**

991126	Microscope Objective Step Down Adapter. Steps down K1 CentriMax S Front \$ for use with Microscope Objective Adapters.	140.00
860600	Z Front for Zeiss ICS objectives. Replaces K 1 CentriMax standard S Front. #991125 is NOT required.	495.00
991125	RMS Objective Adapter (use infinity-corrected type).	150.00
991130	Achrovid Objective Adapter (M26 thread—also suitable for Mitutoyo)	150.00
991131	Nikon CFI60 Objective Adapter (M25 thread)	150.00
860190	M49 Filter Holder.	65.00
	Permits use of M49 photo-type filters on K1 CentriMax MX-1 to MX-5 Objectives.	
770170	Large Mounting Clamp	85.00
770171	Mounting Bar	85.00
770175	ROBUSTO™ Mounting Bar (required for ROBUSTO™)	100.00

Prices and specifications are FOB Centennial, CO and are subject to change without notice.

Model K1 CentriMax, CentriTel, Achrovid, and Nelsonian are trademarks of INFINITY PHOTO-OPTICAL COMPANY.

U.S. Patent 5,452,133. Additional U.S. and foreign patents pending and may apply.